

SACCHAROMYCES Pricing

		Number of Barrels															
		1	2	3	4	5	6	7	8	9	10	15	20	30	40	50	60
Million cells per mL	7	\$104	\$122	\$139	\$162	\$181	\$222	\$239	\$269	\$303	\$336	\$424	\$539	\$662	\$883	\$1,104	1,324
	10	\$111	\$130	\$162	\$207	\$240	\$288	\$336	\$352	\$364	\$404	\$536	\$631	\$946	\$1,261	⊗	⊗
	12	\$116	\$143	\$186	\$231	\$288	\$346	\$355	\$388	\$436	\$485	\$603	\$757	\$1,135	⊗	⊗	⊗
	15	\$130	\$162	\$216	\$288	\$360	\$364	\$424	\$485	\$511	\$536	\$710	\$946	⊗	⊗	⊗	⊗
	20	\$141	\$207	\$288	\$367	\$404	\$485	\$539	\$572	\$589	\$631	\$946	\$1,261	⊗	⊗	⊗	⊗
	25	\$151	\$240	\$360	\$404	\$505	\$536	\$606	\$631	\$710	\$788	\$1,182	⊗	⊗	⊗	⊗	⊗
	30	\$160	\$288	\$372	\$485	\$536	\$612	\$662	\$757	\$851	\$946	\$1,419	⊗	⊗	⊗	⊗	⊗

⊗ Please contact us for a quote. Prices subject to change without notice.

SACCHAROMYCES Pitch Rates

Ales

Although there are many schools of thought on pitch rates, generally a pitch rate of 7 million cells/mL will be adequate for ales under 1.060 OG (15°P). For ales over 15 degrees Plato, we generally recommend a pitch rate of 0.75 mil cells/mL/degree Plato, but this may vary by strain.

Lagers

Recommended pitch rates for lagers depend on pitch temperatures. Ale pitch rates (see above) can be used if pitched in the mid-60s°F (17-19°C) and brought down to desired lager fermentation temperature after activity is observed (typically within 18-24 hours of pitching). If pitching cold around 50°F (10°C), we recommend pitching at a higher rate -- 1.5 million cells/mL/degrees Plato is a good rule of thumb.

Please [contact us](#) with questions or for personalized recommendations.

BRETTANOMYCES Pitch Rates & Pricing

We recommend a pitch rate 500,000 cells/mL for a secondary pitch of Brett.

		Number of Barrels															
		1	2	3	4	5	6	7	8	9	10	15	20	30	40	50	60
Million cells per mL	500k	\$50	\$52	\$61	\$72	\$78	\$84	\$91	\$98	\$106	\$114	\$159	\$167	\$195	\$212	\$226	\$240

LACTOBACILLUS Pitch Rates & Pricing

We recommend direct pitching 1L of Lacto Blend to 3bbls of wort for most kettle sours. If pitched at 95°F (35°C), pH will drop to roughly 3.2-3.4 within 18-48 hours.

Do not let temperature rise above 100°F (38°C) and use only unhopped wort. This Lacto blend is extremely sensitive to hops. Even 1-2 IBUs can inhibit souring.

All bacterial cultures, including the Lacto Blend, are \$82 per liter.